

Oven/Incubator (dual-purpose)

Features

- Polished stainless-steel chamber, semicircular arcs at corners for easy cleaning, and the space between the shelves in the chamber is adjustable.
- Microprocessor temperature controller ensures a precise and reliable control.
- Forced-air convection.
- Double Layer glass door, larger viewing window.
- Independent temperature-limiting Alarm.

Options

- Independent temperature-limiting Alarm system.
- RS485 connector.
- Intelligent programmable temperature controller.



Specifications

Model	PH-010 (A)	PH-030 (A)	PH-050 (A)	PH-070 (A)	PH-140 (A)	PH-240 (A)
Electrical Requirement	220V 50HZ					
Temperature Range	a. as incubator:RT+5~80°C b.as Oven:80~220°C					
Display Resolution	0.1°C					
Temperature Accuracy	a. ±0.5°C b. ±1°C					
Ambient Temperature	+5~35°C					
Power Consumption	850W	850W	1100W	1550W	2050W	2450W
Interior Dimension (W×D×H) mm	250×260×250	340×325×320	420×395×350	450×400×450	550×450×550	600×500×750
Exterior Dimension (W×D×H)mm	530×480×420	620×490×490	720×590×520	740×590×630	840×570×730	880×630×930
Shelves	2(pcs)					
Timing Range	0~5999min					

Hot-air Sterilizer

Microprocessor controller (with timing function)

Summary

- The sterilizer takes advantage of high temperature disinfection principle, which destroys cell bioplasm, denatures protein and concentrates electrolyte to kill microorganisms.

Features

- Microprocessor temperature controller.
- Rapid sterilization time.
- Forced air convection.
- Auto stop when temperature over limit.
- Air inlet adjustment & timing control.



Specifications

Model	GRX-9013A	GRX-9023A	GRX-9053A	GRX-9073A	GRX-9123A	GRX-9203A
Electrical Requirement	220V 50HZ					
Temperature Range	RT+10~200°C / RT+10~250°C					
Temperature Uniformity	0.1°C					
Display Resolution	±1°C					
Temperature Stability	+5~40°C					
Power Consumption	850W	850W	1100W	1550W	2050W	2450W
Interior Dimension (W×D×H) mm	250×260×250	340×325×320	420×395×350	450×400×450	550×450×550	600×550×600
Exterior Dimension (W×D×H)mm	530×480×420	620×490×490	720×590×520	740×590×630	840×570×730	880×750×800
Shelves	2(pcs)					
Timing Range	0~5999min					