

# Two & Four Chambers Heating Incubator

## Features

- Two or four chambers are available, effectively saving laboratory space
- Each chamber is independently controlled, without mutual interference, precise and stable temperature control;
- 304 stainless steel inside , easy to clean and maintain
- Natural heating convection ,no noise, avoiding evaporation of samples during the test.
- With inner glass door for easy observation
- PID controller, accurate temperature control with small fluctuation, with timer function of time 99 hours and 59 minutes.
- Independent over-temperature alarm system ensures experiments running safely.(option)

## Specifications

Model	DHP-150-2B	DHP-100-4B
voltage	AC220V 50HZ	
Temp Range	RT+5~65°C	
Disply Resolution	0.1°C	
Temp Stability	±0.5°C	
Ambient Temperature	+5~35°C	
Power Consumption	1200W	1500W
Chamber Volume	150Lx2	100Lx4
Interior Dimension (WxDxH)mm	560x530x600 /Chamber	455x430x510 /Chamber
External Dimension (WxDxH)mm	705x740x1800	1220x616x1620
Shelves	2PCS (each chamber)	

